





FingerTPS™ Sizing Chart

- 1 Print this page at 100% (do not scale or fit to page). Check scaling with ruler.
- 2 Determine hand size by measuring knuckle width on hand diagram at right.
- 3 Use recommended sensor sizes based on Hand Size Table or match finger/thumb to sensor sizes using the Finger/Thumb Width Chart below.



Hand Size Table

Hand Size	Thumb	Index	Middle	Ring	Pinky
S	TS	S	M	S	XS
M	TM	M	L	M	S
L	TL	L	L	L	M



Align edge of index finger to red line

